

# Spiderflash Recorder Instructions

**Once you are hooked up at the hospital you will never be turning the recorder off**

**DO NOT remove back or change batteries at any time on the spiderflash recorder**

## **How should I wear the Spiderflash?**

- \* Wear the spiderflash recorder in the way most convenient for you – just make sure the green record button is accessible
- \* The spiderflash may be clipped to a belt, slipped in to a shirt pocket

## **Things you can NOT do with the Spiderflash recorder:**

- \* Recorder can not get wet or be exposed to moisture (must remove recorder before bath/shower or swimming)

## **Replacing electrodes:**

- \* Change electrodes every 24-48hrs and every time you bath, shower or swim.
- \* When replacing electrodes make sure skin is washed and dry.
- \* It's okay to move electrode placement around as to help make the skin less irritated
- \* Replace any loose electrodes immediately
- \* Do not use electrodes with dry gel

## **Spiderflash recorder “beep” noises:**

- \* One single beep will occur when you press the green recorder button on the front of the device. This will indicate that the recorder is now saving your recording. It will save a minute and a half before and after you press the button

## **LED flashing light on front of Spiderflash:**

- \* On the front of the spiderflash expect that a green light will flash periodically



Non-Invasive Auto Trigger Event Recorder

SPIDERFLASH-t AFIB

